



## Stock Ampoules


Open and Closed injection ampoules manufactured by SCHOTT Pharmaceutical Packaging from SCHOTT Fliolax® Neutral Type I glass, with Colourbreak or OPC (One Point Cut) for easy opening. Ampoules can be printed and manufactured to your precise specification for viable production runs.

### OPEN AMPOULES - FORM B - FLAME CUT

Ref	Size	Glass Colour	Opening Method	a	b	d	e	n	Pack Size
AIOC001	1ml	Clear	CB	10.75	21	60	8	25.5	646
AIOC002	2ml	Clear	CB	10.75	33	72	8	37.5	646
AIOC005	5ml	Clear	CB	14.75	41	83	9	37.5	343
AIOC010	10ml	Clear	CB	17.75	55	102	9.5	62	240
AIOC020	20ml	Clear	CB	22.5	65	113	11	76	144
AIOA001	1ml	Amber	CB	10.75	21	60	8	25.5	646
AIOA002	2ml	Amber	CB	10.75	33	72	8	37.5	646
AIOA005	5ml	Amber	CB	14.75	41	83	9	37.5	343
AIOA010	10ml	Amber	CB	17.75	55	102	9.5	62	240

### CLOSED AMPOULES - FORM D

Ref	Size	Glass Colour	Opening Method	a	b	d	e	n	Pack Size
AICC001	1ml	Clear	CB	10.75	21	70	6	25.5	646
AICC002	2ml	Clear	CB	10.75	33	83	6	37.5	646
AICC005	5ml	Clear	CB	14.75	41	95	7	46.5	343
AICC010	10ml	Clear	CB	17.75	55	112	7.1	62	240
AICC020	20ml	Clear	CB	22.5	65	126	7.8	76	144
AICC025	25ml	Clear	CB	27.5	107	141	7.8	116.5	144
AICCF001	1ml	Clear	OPC	10.75	21	70	6	25.5	646
AICCF002	2ml	Clear	OPC	10.75	33	83	6	37.5	646
AICCF010	10ml	Clear	OPC	17.75	55	112	7.1	62	240



#### KEY

#### 'I' in the product reference

Indicates that the ampoule conforms to ISO 9187-1 (Colourbreak) ISO 9187-2 (OPC) for injection ampoules

#### Opening Method

CB = Colourbreak - White Breaking  
OPC = One Point Cut - Red Dot

